

USP20 Antibody (R34186)


Catalog No.	Formulation	Size
R34186-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10868
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This USP20 antibody is available for research use only.



IHC testing of FFPE human lymph node with USP20 antibody at 3ug/ml. Microwaved antigen retrieval with Tris/EDTA buffer pH9, HRP-staining.



Western blot testing of human tonsil lysate with USP20 antibody at 1ug/ml. Predicted molecular weight ~102 kDa.

Description

Additional name(s) for this target protein: Ubiquitin specific protease 20

Application Notes

Optimal dilution of the USP20 antibody should be determined by the researcher.

Immunogen

Amino acids HGEQKIEAETRAV were used as the immunogen for this USP20 antibody.

Storage

Aliquot and store the USP20 antibody at -20°C.