


## USP28 Antibody (R36098)

Catalog No.	Formulation	Size
R36098-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57646
Applications	Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This USP28 antibody is available for research use only.



IHC testing of FFPE human testis with USP28 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### Description

Additional name(s) for this target protein: Ubiquitin specific peptidase 28

### Application Notes

Optimal dilution of the USP28 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SVELKHYIQEDN were used as the immunogen for this USP28 antibody.

## **Storage**

Aliquot and store the USP28 antibody at -20oC.