


USP4 Antibody (R33117)

Catalog No.	Formulation	Size
R33117-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7375
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This USP4 antibody is available for research use only.



Western blot testing of human heart lysate with USP4 antibody at 1ug/ml. Predicted molecular weight: 104kDa.

Description

Additional name(s) for this target protein: Ubiquitin specific protease 4

Application Notes

Optimal dilution of the USP4 antibody should be determined by the researcher.

Immunogen

Amino acids AEGGGCRERPDAE were used as the immunogen for this USP4 antibody.

Storage

Aliquot and store the USP4 antibody at -20oC.