


USP6 Antibody (R34484)

Catalog No.	Formulation	Size
R34484-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9098
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This USP6 antibody is available for research use only.



Western blot testing of human placenta lysate with USP6 antibody at 2ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted/observed molecular weight: ~90kDa.

Description

Additional name(s) for this target protein: Ubiquitin specific protease 6, UBI, TRE2, TRE17

Application Notes

Optimal dilution of the USP6 antibody should be determined by the researcher.

Immunogen

Amino acids KISPLHHLQMECSP were used as the immunogen for this USP6 antibody.

Storage

Aliquot and store the USP6 antibody at -20oC.