


## Uteroglobin Antibody / SCGB1A1 [clone AEDA-19] (RQ5374)

Catalog No.	Formulation	Size
RQ5374	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

<b>Availability</b>	1-2 weeks
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	AEDA-19
<b>Purity</b>	Affinity purified
<b>UniProt</b>	P11684
<b>Applications</b>	Western Blot : 1:500-1:2000 Flow Cytometry : 1:50
<b>Limitations</b>	This Uteroglobin antibody is available for research use only.



Western blot testing of human A549 cell lysate with Uteroglobin antibody. Expected molecular weight: 10-16 kDa.

### Description

The SCGB1A1 gene encodes a member of the secretoglobin family of small secreted proteins. Uteroglobin has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma. [RefSeq]

## **Application Notes**

Optimal dilution of the Uteroglobin antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide specific to human Uteroglobin / SCGB1A1 was used as the immunogen for the Uteroglobin antibody.

## **Storage**

Store the Uteroglobin antibody at -20oC.