


Utf1 Antibody (F54761)

Catalog No.	Formulation	Size
F54761-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54761-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity purified
UniProt	Q6J1H4
Applications	Western Blot : 1:500-1:2000
Limitations	This Utf1 antibody is available for research use only.



Western blot testing of mouse 1) brain, 2) liver, 3) heart and 4) testis tissue lysate with Utf1 antibody. Predicted molecular weight: 36 kDa, commonly observed molecular weight: 36-48 kDa.

Description

Acts as a transcriptional coactivator of ATF2. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the Utf1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 325-359 from the mouse protein was used as the immunogen for the Utf1 antibody.

Storage

Aliquot the Utf1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.