


## UXT Antibody Biotin Conjugate (R33472BTN)

Catalog No.	Formulation	Size
R33472BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8409
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:64000
Limitations	This UXT antibody is available for research use only.



Western blot testing of human kidney lysate with biotinylated UXT antibody at 0.3ug/ml.  
Predicted molecular weight: ~20kDa.

### Description

Additional name(s) for this target protein: Ubiquitously-expressed transcript

### Application Notes

Optimal dilution of the UXT antibody should be determined by the researcher.

## **Immunogen**

Amino acids DHRDKVYEQLAKY were used as the immunogen for this UXT antibody.

## **Storage**

Aliquot and store the UXT antibody at -20oC.