


UXT Antibody (R33472)

Catalog No.	Formulation	Size
R33472-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8409
Applications	Western Blot : 0.01-0.03ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This UXT antibody is available for research use only.

kDa 200 –	Western blot testing of human kidney lysate with UXT antibody at 0.03ug/ml. Predicted
116 – 97 –	molecular weight: ~20kDa.
66 –	
44 –	
31 –	
22 –	
14 –	
6-	



IHC staining of FFPE human thyroid gland with UXT antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Cytoplasmic staining with minor nuclear staining is seen.

Description

Additional name(s) for this target protein: Ubiquitously-expressed transcript

Application Notes

Optimal dilution of the UXT antibody should be determined by the researcher.

Immunogen

Amino acids DHRDKVYEQLAKY were used as the immunogen for this UXT antibody.

Storage

Aliquot and store the UXT antibody at -20oC.