

VANGL2 Antibody (R35343)

Catalog No.	Formulation	Size
R35343-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Cow, Human, Zebrafish
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57216
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This VANGL2 antibody is available for research use only.



Western blot testing of mouse brain lysate with VANGL2 antibody at 0.3ug/ml. Predicted molecular weight ~60 kDa.

Description

Additional name(s) for this target protein: Vang-like 2

Application Notes

Optimal dilution of the VANGL2 antibody should be determined by the researcher.

Immunogen

Amino acids QKYLRTTKQQPYH were used as the immunogen for this VANGL2 antibody.

Storage

Aliquot and store the VANGL2 antibody at -20oC.