

VAV3 Antibody (R33336)

Catalog No.	Formulation	Size
R33336-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10451
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This VAV3 antibody is available for research use only.



Western blot testing of human bone marrow lysate with VAV3 antibody at 1ug/ml.
Predicted molecular weight: ~97 kDa (isoform alpha), ~87 kDa (isoform beta).

Description

Additional name(s) for this target protein: Vav 3 guanine nucleotide exchange factor

Application Notes

Optimal dilution of the VAV3 antibody should be determined by the researcher.

Immunogen

Amino acids CSGEQGTKLPEK were used as the immunogen for this VAV3 antibody.

Storage

Aliquot and store the VAV3 antibody at -20oC.