


VCP Antibody (R35027)

Catalog No.	Formulation	Size
R35027-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat, S. cerevisiae
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	851431
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This VCP antibody is available for research use only.



Western blot testing of HeLa lysate with VCP antibody at 0.3ug/ml. Predicted/observed molecular weight: ~89/97kDa.

Description

Additional name(s) for this target protein: Valosin containing protein, CDC48, CDC48p, YDL126C, p97

Application Notes

Optimal dilution of the VCP antibody should be determined by the researcher.

1. This antibody may cross react in many other species.

Immunogen

Amino acids PIADTIEGITGNLFD were used as the immunogen for this VCP antibody.

Storage

Aliquot and store the VCP antibody at -20oC.