


VDAC2 Antibody (R33620)

Catalog No.	Formulation	Size
R33620-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7417
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This VDAC2 antibody is available for research use only.



Western blot testing of human heart lysate with VDAC2 antibody at 0.3ug/ml. Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: Voltage-dependent anion channel 2

Application Notes

Optimal dilution of the VDAC2 antibody should be determined by the researcher.

Immunogen

Amino acids DSAKSKLTRNN were used as the immunogen for this VDAC2 antibody.

Storage

Aliquot and store the VDAC2 antibody at -20oC.