


## VDAC3 Antibody (RQ5850)


| Catalog No. | Formulation                                           | Size   |
|-------------|-------------------------------------------------------|--------|
| RQ5850      | 0.5mg/ml if reconstituted with 0.2ml sterile DI water | 100 ug |

**Bulk quote request**


|                           |                                                                                                                                     |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                                                                                                   |
| <b>Species Reactivity</b> | Human, Mouse, Rat                                                                                                                   |
| <b>Format</b>             | Antigen affinity purified                                                                                                           |
| <b>Host</b>               | Rabbit                                                                                                                              |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                                                                                                          |
| <b>Isotype</b>            | Rabbit IgG                                                                                                                          |
| <b>Purity</b>             | Affinity purified                                                                                                                   |
| <b>Buffer</b>             | Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide                                                                   |
| <b>UniProt</b>            | Q9Y277                                                                                                                              |
| <b>Localization</b>       | Cytoplasmic                                                                                                                         |
| <b>Applications</b>       | Western Blot : 0.5-1ug/ml<br>Immunohistochemistry : 1-2ug/ml<br>Flow Cytometry : 1-3ug/million cells<br>Direct ELISA : 0.1-0.5ug/ml |
| <b>Limitations</b>        | This VDAC3 antibody is available for research use only.                                                                             |




IHC staining of FFPE human breast cancer with VDAC3 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE human testis cancer with VDAC3 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE human melanoma with VDAC3 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE mouse testis with VDAC3 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE rat testis with VDAC3 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) human K562, 2) human A431, 3) rat heart, 4) rat kidney and 5) mouse heart lysate with VDAC3 antibody. Predicted molecular weight: 31 kDa.



Flow cytometry testing of human U-2 OS cells with VDAC3 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= VDAC3 antibody.

## Description

Voltage-dependent anion-selective channel protein 3 (VDAC3) is a protein that in humans is encoded by the VDAC3 gene on chromosome 8. This gene encodes a voltage-dependent anion channel (VDAC), and belongs to the mitochondrial porin family. VDACS are small, integral membrane proteins that traverse the outer mitochondrial membrane and conduct ATP and other small metabolites. They are known to bind several kinases of intermediary metabolism, thought to be involved in translocation of adenine nucleotides, and are hypothesized to form part of the mitochondrial permeability transition pore, which results in the release of cytochrome c at the onset of apoptotic cell death. Alternatively transcript variants encoding different isoforms have been described for this gene.

## Application Notes

Optimal dilution of the VDAC3 antibody should be determined by the researcher.

## Immunogen

Recombinant human protein (amino acids T83-E147) was used as the immunogen for the VDAC3 antibody.

## Storage

After reconstitution, the VDAC3 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.