

VDP p115 Antibody / Vesicle docking protein p115 / USO1 (R33426)

Catalog No.	Formulation	Size
R33426-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8615
Applications	Western Blot:0.3-1ug/ml : ELISA (peptide) LOD : 1:2000
Limitations	This VDP p115 antibody is available for research use only.



Western blot testing of K562 cell lysate with VDP p115 antibody at 0.3ug/ml. Predicted molecular weight: ~108 kDa but routinely observed at 108~115 kDa.

Description

Additional name(s) for this target protein: Vesicle docking protein p115; USO1; General vesicular transport factor p115

Application Notes

Optimal dilution of the VDP p115 antibody should be determined by the researcher.

Immunogen

Amino acids DESEDPGKLDH were used as the immunogen for this VDP p115 antibody.

Storage

Aliquot and store the VDP p115 antibody at -20oC.