


VEGFA Antibody (R35637)

Catalog No.	Formulation	Size
R35637-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7422
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This VEGFA antibody is available for research use only.



Western blot testing of HEK293 lysate with VEGFA antibody at 0.1ug/ml. Expected molecular weight 19~22 kDa (monomer) and 38~44 kDa (dimer). Banding may appear larger than predicted due to glycosylation.

Description

Additional name(s) for this target protein: Vascular endothelial growth factor A

Application Notes

Optimal dilution of the VEGFA antibody should be determined by the researcher.

Immunogen

Amino acids RIKPHQQQHIGE were used as the immunogen for this VEGFA antibody.

Storage

Aliquot and store the VEGFA antibody at -20oC.