


VHL Antibody Biotin Conjugate (R35773BTN)

Catalog No.	Formulation	Size
R35773BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7428
Applications	Western Blot : 0.3-1ug/ml Enzyme ImmunoAssay : suitable ELISA (peptide) LOD : 1:128000
Limitations	This VHL antibody is available for research use only.



Western blot testing of human ovary lysate with biotinylated VHL antibody at 0.3ug/ml.
Predicted molecular weight ~23/19kDa (Isoforms 1/2).

Description

Additional name(s) for this target protein: Von Hippel-Lindau tumor suppressor, VHL1

Application Notes

Optimal dilution of the VHL antibody should be determined by the researcher.

1. Isoforms 1 (213 AA) and 2 (172 AA) may be detected with this Ab.

Immunogen

Amino acids RSLVKPENYRRLD were used as the immunogen for this VHL antibody.

Storage

Aliquot and store the VHL antibody at -20oC.