


Vigilin Antibody / HDLBP / High density lipoprotein-binding protein [clone 28H81] (RQ8530)

Catalog No.	Formulation	Size
RQ8530	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28H81
Purity	Affinity chromatography
UniProt	Q00341
Applications	Western Blot : 1:500-1:2000
Limitations	This Vigilin antibody is available for research use only.



Western blot testing of human HeLa cell lysate with Vigilin antibody. Predicted molecular weight ~141 kDa.

Description

Appears to play a role in cell sterol metabolism. It may function to protect cells from over-accumulation of cholesterol.
[UniProt]

Application Notes

Optimal dilution of the Vigilin antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human Vigilin protein was used as the immunogen for the Vigilin antibody.

Storage

Store the Vigilin antibody at -20oC.