


VLCAD Antibody (R33512)

Catalog No.	Formulation	Size
R33512-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	37
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This VLCAD antibody is available for research use only.



Western blot testing of human heart lysate with VLCAD antibody at 2ug/ml. Predicted molecular weight: ~68kDa.

Description

Additional name(s) for this target protein: Acyl-Coenzyme A dehydrogenase, very long chain; ACADVL

Application Notes

Optimal dilution of the VLCAD antibody should be determined by the researcher.

Immunogen

Amino acids DKSDSHPSDALTRK were used as the immunogen for this VLCAD antibody.

Storage

Aliquot and store the VLCAD antibody at -20oC.