


## VPS29 Antibody (R33404)

Catalog No.	Formulation	Size
R33404-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51699
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This VPS29 antibody is available for research use only.



Western blot testing of human spleen lysate with VPS29 antibody at 0.1ug/ml. Predicted molecular weight: ~20/21kDa (isoforms 1/2).

## Description

Additional name(s) for this target protein: Vacuolar protein sorting 29

## Application Notes

Optimal dilution of the VPS29 antibody should be determined by the researcher.

## Immunogen

Amino acids DDVKVERIEYKKP were used as the immunogen for this VPS29 antibody.

## Storage

Aliquot and store the VPS29 antibody at -20oC.