


VPS37C Antibody (R34021)


| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34021-100UG | 0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|--|
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 55048 |
| Applications | Western Blot: 0.03-0.1ug/ml Immunohistochemistry (FFPE): 3-6ug/ml ELISA (peptide) LOD: 1:32000 |
| Limitations | This VPS37C antibody is available for research use only. |



IHC staining of FFPE human prostate with VPS37C antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



Western blot testing of human duodenum lysate with VPS37C antibody at 0.03g/ml. Predicted/observed molecular weight: ~39kDa.

Description

Additional name(s) for this target protein: Vacuolar protein sorting 37C

Application Notes

Optimal dilution of the VPS37C antibody should be determined by the researcher.

Immunogen

Amino acids ETLKDKTLQELEELQ were used as the immunogen for this VPS37C antibody.

Storage

Aliquot and store the VPS37C antibody at -20oC.