


## VPS37C Antibody (R34021)


Catalog No.	Formulation	Size
R34021-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55048
Applications	Western Blot : 0.03-0.1ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This VPS37C antibody is available for research use only.



IHC staining of FFPE human prostate with VPS37C antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



Western blot testing of human duodenum lysate with VPS37C antibody at 0.03g/ml. Predicted/observed molecular weight: ~39kDa.

## Description

Additional name(s) for this target protein: Vacuolar protein sorting 37C

## Application Notes

Optimal dilution of the VPS37C antibody should be determined by the researcher.

## Immunogen

Amino acids ETLKDKTLQEEELQ were used as the immunogen for this VPS37C antibody.

## Storage

Aliquot and store the VPS37C antibody at -20oC.