


## VPS41 Antibody (R33468)

Catalog No.	Formulation	Size
R33468-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27072
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This VPS41 antibody is available for research use only.



IHC testing of FFPE human kidney with VSP41 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

## Description

Additional name(s) for this target protein: Vacuolar protein sorting 41

## Application Notes

Optimal dilution of the VPS41 antibody should be determined by the researcher.

## Immunogen

Amino acids DHLKKDSQNKTLTK were used as the immunogen for this VPS41 antibody.

## **Storage**

Aliquot and store the VPS41 antibody at -20oC.