


VPS45 Antibody (R34002)

Catalog No.	Formulation	Size
R34002-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11311
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This VPS45 antibody is available for research use only.



Western blot testing of human peripheral blood mononucleocyte lysate with VPS45 antibody at 0.3ug/ml. Predicted molecular weight: ~65kDa.

Description

Additional name(s) for this target protein: Vacuolar protein sorting 45

Application Notes

Optimal dilution of the VPS45 antibody should be determined by the researcher.

Immunogen

Amino acids ESSQVTSRSASRR were used as the immunogen for this VPS45 antibody.

Storage

Aliquot and store the VPS45 antibody at -20°C.