


VPS45 Antibody (R34063)


Catalog No.	Formulation	Size
R34063-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11311
Localization	Cytoplasmic
Applications	Western Blot : 0.03-0.1 ug/ml IHC (FFPE) : 2.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This VPS45 antibody is available for research use only.



Western blot testing of human placenta lysate with VPS45 antibody at 0.03ug/ml.
Predicted molecular weight: ~65 kDa.



IHC testing of FFPE human spleen tissue with VPS45 antibody at 2.5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Vacuolar protein sorting 45

Application Notes

Optimal dilution of the VPS45 antibody should be determined by the researcher.

Immunogen

Amino acids FQKKKPKEQQKLES were used as the immunogen for this VPS45 antibody.

Storage

Aliquot and store the VPS45 antibody at -20oC.