


VRK2 Antibody (R35087)

Catalog No.	Formulation	Size
R35087-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7444
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This VRK2 antibody is available for research use only.



Western blot testing of K562 nuclear lysate with VRK2 antibody at 2ug/ml. Predicted molecular weight: ~58kDa.

Description

Additional name(s) for this target protein: Vaccinia related kinase 2

Application Notes

Optimal dilution of the VRK2 antibody should be determined by the researcher.

Immunogen

Amino acids PNSQKVDSQKAATK were used as the immunogen for this VRK2 antibody.

Storage

Aliquot and store the VRK2 antibody at -20oC.