


VSV-G Tag Antibody [clone DB-22] (RQ4681)

Catalog No.	Formulation	Size
RQ4681	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	VSV-G Tag
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	DB-22
Purity	Affinity purified
Applications	Western Blot : 1:3000-1:10000
Limitations	This VSV-G Tag antibody is available for research use only.



Western blot testing of a VSV-G tagged fusion protein with VSV-G Tag antibody.

Description

Epitope tags consisting of short sequences recognized by well-characterized antibodies have been widely used in the study of protein expression in various systems.

Application Notes

Optimal dilution of the VSV-G Tag antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to VSV-G was used as the immunogen for the VSV-G Tag antibody.

Storage

Store the VSV-G Tag antibody at -20oC.