


## VTCN1 Antibody (R35006)

Catalog No.	Formulation	Size
R35006-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79679
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This VTCN1 antibody is available for research use only.



Western blot testing of human spleen lysate with VTCN1 antibody at 0.2ug/ml. Predicted molecular weight: ~31kDa.

## Description

Additional name(s) for this target protein: V-set domain containing T cell activation inhibitor 1

## Application Notes

Optimal dilution of the VTCN1 antibody should be determined by the researcher.

## Immunogen

Amino acids PEVNVDYNASSET were used as the immunogen for this VTCN1 antibody.

## **Storage**

Aliquot and store the VTCN1 antibody at -20oC.