


VTI1B Antibody (R35555)

Catalog No.	Formulation	Size
R35555-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10490
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This VTI1B antibody is available for research use only.



Description

Additional name(s) for this target protein: Vesicle transport through interaction with t-SNAREs 1B

Application Notes

Optimal dilution of the VTI1B antibody should be determined by the researcher.

Immunogen

Amino acids RDFDEKQQEANET were used as the immunogen for this VTI1B antibody. **Storage** Aliquot and store the VTI1B antibody at -20oC.