


VT11B Antibody (R35555)

| Catalog No. | Formulation | Size |
|--------------|-------------------------------------------------------------------------------|--------|
| R35555-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|----------------------------------------------------------|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Dog, Mouse, Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 10490 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This VT11B antibody is available for research use only. |



Western blot testing of human duodenum lysate with VT11B antibody at 1ug/ml.
Predicted molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: Vesicle transport through interaction with t-SNAREs 1B

Application Notes

Optimal dilution of the VT11B antibody should be determined by the researcher.

Immunogen

Amino acids RDFDEKQQEANET were used as the immunogen for this VT11B antibody.

Storage

Aliquot and store the VT11B antibody at -20oC.