


WNT16 Antibody [clone 28W88] (RQ8627)

Catalog No.	Formulation	Size
RQ8627	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28W88
Purity	Affinity chromatography
UniProt	Q9UBV4
Applications	Western Blot : 1:500-1:2000
Limitations	This WNT16 antibody is available for research use only.



Western blot testing of human HepG2 cell lysate with WNT16 antibody. Predicted molecular weight ~41 kDa.

Description

WNT16 protein is a ligand for members of the frizzled family of seven transmembrane receptors. [UniProt]

Application Notes

Optimal dilution of the WNT16 antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human Protein Wnt-16 was used as the immunogen for the WNT16 antibody.

Storage

Store the WNT16 antibody at -20oC.