

# XAB2 Antibody (R36101)


Catalog No.	Formulation	Size
R36101-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	56949
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 3-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This XAB2 antibody is available for research use only.



Western blot testing of A431 lysate with XAB2 antibody at 0.3ug/ml. Predicted molecular weight: ~100kDa.



IHC staining of FFPE human cerebral cortex with XAB2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## **Description**

Additional name(s) for this target protein: XPA binding protein 2, HCNP, HCRN

## **Application Notes**

Optimal dilution of the XAB2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SVPAAVFGSLKED were used as the immunogen for this XAB2 antibody.

## **Storage**

Aliquot and store the XAB2 antibody at -20oC.