


XAGE1 Antibody (isoform d) (R33236)

Catalog No.	Formulation	Size
R33236-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9503
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This XAGE1 antibody is available for research use only.



IHC testing of FFPE human testis with XAGE1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: G antigen, family D, 2; GAGED2

Application Notes

Optimal dilution of the XAGE1 antibody should be determined by the researcher.

1. This XAGE1 antibody is specific for isoform d only.

Immunogen

Amino acids CGFGFRRQGEDNT were used as the immunogen for this XAGE1 antibody.

Storage

Aliquot and store the XAGE1 antibody at -20oC.