


YAP1 Antibody (R36318)

Catalog No.	Formulation	Size
R36318-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10413
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This YAP1 antibody is available for research use only.



Western blot testing of MDCK lysate with YAP1 antibody at 0.5ug/ml. Predicted molecular weight: 54 kDa but routinely observed at 65-70 kDa.

Description

Additional name(s) for this target protein: Yes-associated protein 1; YAP65

Application Notes

Optimal dilution of the YAP1 antibody should be determined by the researcher.

Immunogen

Amino acids TYHSRDESTDS were used as the immunogen for this YAP1 antibody.

Storage

Aliquot and store the YAP1 antibody at -20oC.