


YWHAZ Antibody / 14-3-3 zeta/delta [clone 6G5] (RQ5642)


Catalog No.	Formulation	Size
RQ5642	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Monkey
Format	Antigen affinity purified
Host	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Clone Name	6G5
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	P63104
Applications	Western Blot : 0.5-1ug/ml Flow Cytometry : 1-3ug/million cells
Limitations	This YWHAZ antibody is available for research use only.



Flow cytometry testing of human PC-3 cells with YWHAZ antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= YWHAZ antibody.



Flow cytometry testing of human SiHa cells with YWHAZ antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= YWHAZ antibody.



Western blot testing of 1) human HeLa, 2) human A549, 3) monkey COS-7, 4) human Raji, 5) human Caco-2, 6) human Jurkat, 7) mouse brain and 8) rat brain lysate with YWHAZ antibody. Predicted molecular weight ~28 kDa.

Description

14-3-3 protein zeta/delta is a protein that in humans is encoded by the YWHAZ gene on chromosome 8. This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.

Application Notes

Optimal dilution of the YWHAZ antibody should be determined by the researcher.

Immunogen

Amino acids LLEKFLIPNASQAESKVFYFLKMKGDYRYLAEVAAGDDKKGIVDQ from the human protein were used as the immunogen for the YWHAZ antibody.

Storage

After reconstitution, the YWHAZ antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.