


ZAP70 Antibody [clone RM408] (R20424)

Catalog No.	Formulation	Size
R20424-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM408
Purity	Protein A purified from animal origin-free supernatant
UniProt	P43403
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:100-1:200 Western Blot : 1:250-1:500
Limitations	This recombinant ZAP70 antibody is available for research use only.



Western blot testing of human Jurkat cell lysate with recombinant ZAP70 antibody at 1:500 dilution. Predicted molecular weight ~70 kDa.



IHC staining of FFPE human tonsil tissue with recombinant ZAP70 antibody at 1:200.

Description

This antibody reacts to human ZAP70.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant ZAP70 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to residues near the N-terminus of human ZAP70 was used as the immunogen for the recombinant ZAP70 antibody.

Storage

Store the recombinant ZAP70 antibody at -20oC.