


ZBTB24 Antibody (R35173)

Catalog No.	Formulation	Size
R35173-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9841
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This ZBTB24 antibody is available for research use only.



Western blot testing of Daudi lysate with ZBTB24 antibody at 1ug/ml. Predicted molecular weight: ~78kDa.

Description

Additional name(s) for this target protein: Zinc finger and BTB domain containing 24

Application Notes

Optimal dilution of the ZBTB24 antibody should be determined by the researcher.

Immunogen

Amino acids PVEKDENYDPKTED were used as the immunogen for this ZBTB24 antibody.

Storage

Aliquot and store the ZBTB24 antibody at -20oC.