


Zcchc11 Antibody (R35337)


| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35337-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Mouse, Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 230594 |
| Applications | Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This Zcchc11 antibody is available for research use only. |



Western blot testing of rat brain lysate with Zcchc11 antibody at 0.05ug/ml.
Predicted/observed molecular weight (rat): ~163kDa.



Western blot of different testis lysates from wildtype (wt) 8-9 week old C57BL/6 and knock-out (KO) mice with Zcchc11 antibody at 0.2ug/ml.

Description

Additional name(s) for this target protein: Zinc finger, CCHC domain containing 11, Terminal Uridylyltransferase 4

Application Notes

Optimal dilution of the Zcchc11 antibody should be determined by the researcher.

Immunogen

Amino acids QKSAEKQGDLTPSET were used as the immunogen for this Zcchc11 antibody.

Storage

Aliquot and store the Zcchc11 antibody at -20oC.