


# ZDHHC2 Antibody / REAM / ZNF372 (F54830)


Catalog No.	Formulation	Size
F54830-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54830-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
UniProt	Q9UIJ5
Localization	Cytoplasmic
Applications	Western Blot : 1:500-1:1000 Immunohistochemistry (FFPE) : 1:50-1:100
Limitations	This ZDHHC2 antibody is available for research use only.



Western blot testing of human Ramos cell lysate in the absence and presence of immunizing/blocking peptide with ZDHHC2 antibody. Predicted molecular weight ~42 kDa



IHC testing of FFPE human brain tissue with ZDHHC2 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.

## **Description**

Palmitoyltransferase that catalyzes the addition of palmitate onto various protein substrates and is involved in a variety of cellular processes. [UniProt]

#### **Application Notes**

The stated application concentrations are suggested starting points. Titration of the ZDHHC2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

#### **Immunogen**

A portion of amino acids 77-106 from the human protein was used as the immunogen for the ZDHHC2 antibody.

### **Storage**

Aliquot the ZDHHC2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.