


ZDHHC2 Antibody / REAM / ZNF372 (F54830)


Catalog No.	Formulation	Size
F54830-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54830-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
UniProt	Q9UIJ5
Localization	Cytoplasmic
Applications	Western Blot : 1:500-1:1000 Immunohistochemistry (FFPE) : 1:50-1:100
Limitations	This ZDHHC2 antibody is available for research use only.



Western blot testing of human Ramos cell lysate in the absence and presence of immunizing/blocking peptide with ZDHHC2 antibody. Predicted molecular weight ~42 kDa.



IHC testing of FFPE human brain tissue with ZDHHC2 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.

Description

Palmitoyltransferase that catalyzes the addition of palmitate onto various protein substrates and is involved in a variety of cellular processes. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the ZDHHC2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 77-106 from the human protein was used as the immunogen for the ZDHHC2 antibody.

Storage

Aliquot the ZDHHC2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.