

Zebrafish Patched Antibody (R35480)

Catalog No.	Formulation	Size
R35480-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q98864
Applications	Immunofluorescence : suitable (1) ELISA (peptide) LOD : 1:32000
Limitations	This Zebrafish Patched antibody is available for research use only.



Description

Zebrafish Patched Antibody detects Patched protein.

This Zebrafish antibody is part of a [broader Zebrafish / Danio rerio antibody panel](#) offered by NSJ Bioreagents.

Application Notes

Optimal dilution of the Zebrafish Patched antibody should be determined by the researcher.

1. Has been used to study tumors in zebrafish (PMID: 19555497)

Immunogen

Amino acids QTGSKKEPFNYSQ were used as the immunogen for this Zebrafish Patched antibody.

Storage

Aliquot and store the Zebrafish Patched antibody at -20oC.