


## Zebrafish Thrβ Antibody (R33481)

Catalog No.	Formulation	Size
R33481-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Zebrafish
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Applications</b>	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This zebrafish Thrβ antibody is available for research use only.



Western blot testing of three different zebrafish adult heart lysate samples with Thrβ antibody at 2.5ug/ml. Predicted molecular weight: ~44kDa.

### Description

Zebrafish Thrβ Antibody detects Thyroid hormone receptor beta protein.

This Zebrafish antibody is part of a [broader Zebrafish / Danio rerio antibody panel](#) offered by NSJ Bioreagents.

### Application Notes

Optimal dilution of the Zebrafish Thrβ antibody should be determined by the researcher.

## **Immunogen**

Amino acids DIGSAPIVNAPEGNKVD were used as the immunogen for this zebrafish Thrb antibody.

## **Storage**

Aliquot and store the Zebrafish Thrb antibody at -20oC.