


Zebrafish Thrβ Antibody (R33481)

Catalog No.	Formulation	Size
R33481-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This zebrafish Thrβ antibody is available for research use only.

50 kd
37 kd



Western blot testing of three different zebrafish adult heart lysate samples with Thrβ antibody at 2.5ug/ml. Predicted molecular weight: ~44kDa.

Description

Additional name(s) for this target protein: Thyroid hormone receptor beta

Application Notes

Optimal dilution of the Thrβ antibody should be determined by the researcher.

Immunogen

Amino acids DIGSAPIVNAPEGNKVD were used as the immunogen for this zebrafish Thrβ antibody.

Storage

Aliquot and store the Thrb antibody at -20oC.