


Zfp157 Antibody (R35208)

Catalog No.	Formulation	Size
R35208-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	O88219
Gene ID	72154
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Zfp157 antibody is available for research use only.



Western blot testing of mouse spinal cord lysate with Zfp157 antibody at 1ug/ml.
Predicted molecular weight: ~66kDa.

Description

Additional name(s) for this target protein: Zinc finger protein 157; Mszf12

Application Notes

Optimal dilution of the Zfp157 antibody should be determined by the researcher.

Immunogen

Amino acids QCSKQHAYYPSVKP were used as the immunogen for this Zfp157 antibody.

Storage

Aliquot and store the Zfp157 antibody at -20oC.