


ZFP91 Antibody / Zinc finger protein 91 (F55085)

Catalog No.	Formulation	Size
F55085-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F55085-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-2 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
UniProt	Q96JP5
Applications	Western Blot : 1:1000-1:2000
Limitations	This ZFP91 antibody is available for research use only.



Western blot testing of human HL60 cell lysate with ZFP91 antibody. Predicted molecular weight ~63 kDa and ~60 kDa (two isoforms).



Western blot testing of human testis tissue lysate with ZFP91 antibody. Predicted molecular weight ~63 kDa and ~60 kDa (two isoforms).

Description

ZFP91 belongs to a family of proteins called zinc finger proteins, which are known for their ability to bind to specific DNA sequences and regulate gene expression. ZFP91 contains zinc finger domains that enable it to target and bind to DNA, allowing it to modulate the activity of genes involved in biological processes. Studies have shown that ZFP91 plays a crucial role in diverse cellular functions, such as cell growth, proliferation, and differentiation. Additionally, research has revealed that ZFP91 is involved in the regulation of immune responses and inflammation, highlighting its significance in maintaining a healthy immune system.

Application Notes

Titration of the ZFP91 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 216-245 from the human protein was used as the immunogen for this ZFP91 antibody.

Storage

Aliquot the ZFP91 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.