


ZHX2 Antibody / Zinc fingers and homeoboxes protein 2 [clone 28Z52] (RQ8400)

Catalog No.	Formulation	Size
RQ8400	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28Z52
Purity	Affinity chromatography
UniProt	Q9Y6X8
Applications	Western Blot : 1:500-1:2000
Limitations	This ZHX2 antibody is available for research use only.



Western blot testing of human HeLa cell lysate with ZHX2 antibody. Predicted molecular weight ~92 kDa, commonly observed at 92-110 kDa.

Description

ZHX2 acts as a transcriptional repressor, repressing the promoter activity of the CDC25C gene stimulated by NFYA. [UniProt]

Application Notes

Optimal dilution of the ZHX2 antibody should be determined by the researcher.

Immunogen

A synthetic peptide derived from human ZHX2 was used as the immunogen for the ZHX2 antibody.

Storage

Store the ZHX2 antibody at -20oC.