


ZIC2 Antibody (R34286)

Catalog No.	Formulation	Size
R34286-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7546
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This ZIC2 antibody is available for research use only.



Description

Additional name(s) for this target protein: Zic family member 2

Application Notes

Optimal dilution of the ZIC2 antibody should be determined by the researcher.

Immunogen

Amino acids HSGLSSNFNEWY were used as the immunogen for this ZIC2 antibody.

Storage

Aliquot and store the ZIC2 antibody at -20oC.