


ZNF217 Antibody (R34536)

Catalog No.	Formulation	Size
R34536-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7764
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ZNF217 antibody is available for research use only.



IHC testing of FFPE human prostate with ZNF217 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Nuclear staining in select secretory cells is seen.

Description

Additional name(s) for this target protein: Zinc finger protein 217, ZABC1

Application Notes

Optimal dilution of the ZNF217 antibody should be determined by the researcher.

Immunogen

Amino acids KNAADDSDADKVNKN were used as the immunogen for this ZNF217 antibody.

Storage

Aliquot and store the ZNF217 antibody at -20oC.