


## ZNF217 Antibody (R34536)

Catalog No.	Formulation	Size
R34536-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7764
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ZNF217 antibody is available for research use only.



IHC testing of FFPE human prostate with ZNF217 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Nuclear staining in select secretory cells is seen.

### Description

Additional name(s) for this target protein: Zinc finger protein 217, ZABC1

### Application Notes

Optimal dilution of the ZNF217 antibody should be determined by the researcher.

### Immunogen

Amino acids KNAADDSDADKVNKN were used as the immunogen for this ZNF217 antibody.

## Storage

Aliquot and store the ZNF217 antibody at -20oC.