


## ZNF652 Antibody (R35628)

Catalog No.	Formulation	Size
R35628-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22834
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ZNF652 antibody is available for research use only.



Western blot testing of human heart lysate with ZNF652 antibody at 0.01ug/ml. Predicted molecular weight: ~70kDa, routinely observed at 85-90kDa. (Ref 1)

## Description

Additional name(s) for this target protein: Zinc finger protein 652

## Application Notes

Optimal dilution of the ZNF652 antibody should be determined by the researcher.

## Immunogen

Amino acids RHLEKNSSAQHH were used as the immunogen for this ZNF652 antibody.

## Storage

Aliquot and store the ZNF652 antibody at -20oC.

## References (1)