


ZRANB1 Antibody (R34697)

Catalog No.	Formulation	Size
R34697-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54764
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This ZRANB1 antibody is available for research use only.



Western blot testing of human cerebellum lysate with ZRANB1 antibody at 0.1 ug/ml.
Predicted/observed molecular weight: ~80kDa.

Description

Additional name(s) for this target protein: Zinc finger, RAN-binding domain containing 1 protein, TRABID

Application Notes

Optimal dilution of the ZRANB1 antibody should be determined by the researcher.

Immunogen

Amino acids KWLDRYRQIRPCTS were used as the immunogen for this ZRANB1 antibody.

Storage

Aliquot and store the ZRANB1 antibody at -20oC.